

IVD

Fecal Occult Blood Test Cassette

SENSATEST® One Step Fecal Occult Blood Test

“Fecal occult blood” is generally defined as a blood loss of less than 50mL/d. Most of the diseases can cause hidden blood in the stool.

In the early stages, gastrointestinal problems such as colon cancer, ulcers, polyps, colitis, diverticulitis, and fissures may not show any visible symptoms, only occult blood.

Among these diseases, colorectal cancer is the third most common cancer in the world.

Specification

The **SENSATEST®** One Step Fecal Occult Blood Test Device (Faeces) is a rapid test to qualitatively detect lower levels of fecal occult blood in faeces. The test is a qualitative, lateral flow immunoassay for the detection of above 50ng/mL occult blood in Human faeces.

- ◆ Early diagnosis of colon cancer
- ◆ Specimen : Human faeces
- ◆ Test result in 5 minutes


Storage & Stability

Stored at 4-30 °C. Kit is stable up to expiry date if kept in desired storage tempera-

Materials provided

- ◆ Test device
- ◆ Specimen collection tube with extraction buffer
- ◆ Package insert

Procedure



Product Code	Test Medium	Test Format	Packing
FOB 1102	Faeces	Cassette	1 x 40

For further information contact us at :

NEOPHARMA BIOTECH ASIA SDN BHD (648086-W)

Suite 410, Block A, Glomac Business Centre, 10 Jalan SS 6/1, Kelana Jaya, 47301 Petaling Jaya, Selangor Darul Ehsan, Malaysia.

Website : www.neopharmaasia.com Tel : +603-7880 9787 Fax : +603-7880 9296