



## IVD Infectious Disease Tests

### One Step Hepatitis B Surface Antigen Test (Blood/Serum/Plasma)

One Step HBsAg Test is a qualitative immunoassay that detects Hepatitis B Surface Antigen (HBsAg) in serum, plasma or blood samples.


#### SPECIFICATIONS

- Qualitative results within 10-20 minutes.
- 1 ng/mL sensitivity.


#### CERTIFICATIONS

- Institute for Medical Research Malaysia (IMR)
- Certificate of Free Sale (CFS)
- Registered under Medical Devices Establishments (MeDVER)
- CE / ISO 9001

This test is available in strip and cassette format. The procedure is simple and requires no additional instruments to read the results. The **One-Step HBsAg Test** is an easy and effective aid for the screening and diagnosis of Hepatitis B infection.



Cassette Format



Strip Format

### Interpretation of Results

**Negative** : Only one coloured band appears on the control (C) region. No apparent band on the test (T) region. A negative result indicates that HBs Ag is not detected

**Positive** : In addition to a pink coloured control (C) band, a distinct pink coloured band will also appear in the test (T) region. A positive result indicates presence of HBsAg. Please consult your physician to perform a much more detailed examination.

**Invalid** : A total absence of colour in both regions is an indication of procedure error or the test reagent has deteriorated. Repeat the test with a new test device.

**Effectiveness : 99%**

| Product Code | Test Medium            | Test Format | Packing | Store |
|--------------|------------------------|-------------|---------|-------|
| HBsAg 210    | Serum / Plasma / Blood | Strip       | 1 x 100 |       |
| HBsAg 211    | Serum / Plasma / Blood | Cassette    | 1 x 40  |       |

For further information contact us at :

**NEOPHARMA BIOTECH ASIA SDN BHD** (648086-W)

Suite 410, Block A, Glomac Business Centre, 10 Jalan SS 6/1, Kelana Jaya, 47301 Petaling Jaya, Selangor Darul Ehsan, Malaysia.

Website : [www.neopharmaasia.com](http://www.neopharmaasia.com) Tel : +603-7880 9787 Fax : +603-7880 9296