

# Vaperon™ SARS-CoV-2 Neutralizing Antibody Rapid Detection Kit

## ◆ Material Provided:


- Test Cassettes (20 ×)
- Assay Diluent [6ml/bottle] (2 ×)
- Alcohol Prep Pad (20 ×)
- Sterile Disposable Lancet (20 ×)
- Blood Transfer Pipette (20 ×)
- Instruction for Use (1 ×)



## ◆ Intended Use:

Vaperon™ SARS-CoV-2 Neutralizing Antibody Rapid Detection Kit is intended for identifying individuals with adaptive immune response to SARS-CoV-2, indicating recent or prior infection or vaccination.

## ◆ Test Procedures and Result Interpretations:



## ◆ Order Reference Information:

Cat.No	Description of Products	Format	Packing Size	Specimen
COVNA-503	SARS-CoV-2 Neutralizing Antibody Rapid Detection Kit	Cassette	20T	S/P/WB



**Atlas-Link Biotech**

Atlas Link Technology Co., Ltd  
Gu'an South Industry Zone, 065500 Langfang City,  
Hebei Province, the People's Republic of China

Authorised Representative:

**NEOPHARMA BIOTECH ASIA SDN BHD**  
Unit 3.1, Block A, Dataran PHB, Persiaran Golf,  
Saujana Resort, Sek U2, 40150 Shah Alam,  
Selangor Malaysia  
Email: [info@neopharmaasia.com](mailto:info@neopharmaasia.com)